



Monitoring system for Hyperbaric Chambers

GDA monitoring system is available for integration with any make and model of the hyperbaric chambers. This system simplifies and increases the efficiency of chamber operation, and is easy to understand and use.

In order to achieve maximum patient safety. A checklist feature prevents chamber operation unless all start-up procedures are correctly followed. The log book section provides space to enter comments and the system offers a number of different treatment tables to choose from.



The following parameters are recorded:

- Depth
- Time
- Oxygen
- Carbon-dioxide
- Humidity
- Temperature

- Monitoring of what medium that is connected to the bibs.
- Choice between various treatment tables
- Log book for entering comments
- Checklist for patient safety and a fail-safe chamber that will not operate in the event of a check list error

The InData dialog box contains the following fields and controls:

- File name: NEW FILE
- Doctor in charge:
- Nurse:
- Operators:
- Table: TABEL CHOICE
- Buttons: OK, Cancel